

iReach Magnum

Powered Endoscopic
Linear Stapler

55°

55° articulation



Powered
12 Firings



Superior
Stability



Superior
Compression



New blade in
every cartridge

ENDO III Platform powered up

Fundamental essentials of ENDO III stapling system now re-imagined and reassembled into a brand-new entity.

Signature tissue compression paired with smooth powered firing yields unprecedented experience, unlocking new possibilities in the OR.

Superior compression combined with maximum safety



Superior compression that maximize compression, reduce tissue damage and better staple formation




Audio feedback when **15-seconds** precompression is achieved.



Safety release button, reverse switch and manual reverse guarantee the maximum safety.

Built to perform



Up to **55°** of articulation in both directions in 5 gears.



Advanced algorithm on the chip precisely modulate power output to achieve optimal effects



Each battery prepackaged with the handle allows up to 12 times smooth firing.



Indication light shows the status of the battery




Use articulation lever to archive articulation with out the assistance of tissues.

IM45AS/IM45AM/IM45AL compatible with



IM60AS/IM60AM/IM60AL compatible with



ENDO III Cartridge



A wide range of single-patient-use reloads that offers excellent and smooth tissue dividing experience, compatible with 45 or 60 mm staple cartridge.



Cartridge variety
6 staple heights in both 45 and 60 mm length are available for any surgical situation



Single-use blade in every cartridge
Ensures clean cuts and minimizes tumor metastasis.



All cartridges, including green and black through 12 mm trocars

Product Code	Cartridge Color	Open Staple Height	Closed Staple Height	Tissue Thickness	Number of staples	Number of staple rows	QTY/ BOX
REC 45 GRA	Gray	2.0 mm	0.75 mm	Vascular/Thin	66	6	12
REC 60 GRA					90		
REC 45 WHT	White	2.5 mm	1.0 mm	Vascular/Thin	66		
REC 60 WHT					90		
REC 45 BLU	Blue	3.5 mm	1.5 mm	Regular	66		
REC 60 BLU					90		
REC 45 GLD	Gold	3.8 mm	1.75 mm	Regular	66		
REC 60 GLD					90		
REC 45 GRN	Green	4.1 mm	2.0 mm	Regular/Thick	66		
REC 60 GRN					90		
REC 45 BLK	Black	4.4 mm	2.2 mm	Regular/Thick	66		
REC 60 BLK					90		

CONTACT US

202101b

Global Headquarter

120 Xinxing Road
West Zone, TEDA
300462 Tianjin
China

T +86-22-25335069

R&D Center

4480 Lake Forest Dr, Suite 414
Cincinnati
45242 OH
USA

T +1(513)769-1918

F +1(513)769-1921

International Sales Office

Room 1011, 228 Mei Yuan Rd
Enterprise Square, Jing'an District
200040, Shanghai
China

T +86-21-60126560